



ADDITIONAL INFORMATION

Weight	2.48 kg
Dimensions	35 × 21.5 × 43.5 cm
Carton Qty	1
Pricing Unit	Each
Size	41 x 33 x 19cm
Min Order Qty	1
Tax Schedule	Taxable
Item Length	41
Item Width	33
Item Height	19
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The "Heart of the Nation" labeled CardiAct CC-90 Alarmed Outdoor Cabinet provides a compact and secure cabinet for storage of most AEDs. Designed for confined spaces, with its compact size the cabinet also provides clear visibility of the AED with a transparent polycarbonate door that includes AED signage for easy identification in emergencies. Defending against the weather, the CC-90 provides weatherproof protection against dust, rain, cold, rust and humidity with an IP rating of IP54. The CC-90 is suitable for use with the CardiAct Cabinet Stand for versatile mounting options.

The CC-90s Yellow design and Heart of the Nation labelling enhances visibility making it an ideal choice for AEDs located in the public and accessible to anyone. This cabinet ensures that life-saving Automated External Defibrillators (AEDs) are easily accessible, promoting rapid response.

IP Rating: IP54



HOTN CARDIAC ALARMED AED CABINET (YELLOW)

HOTN-CC-90