



## ADDITIONAL INFORMATION

<b>Weight</b>	2 kg
<b>Dimensions</b>	41 × 33 × 19 cm
<b>Pricing Unit</b>	Each
<b>Size</b>	41 x 33 x 19cm
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	31
<b>Item Width</b>	33
<b>Item Height</b>	19
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



## PRODUCT DESCRIPTION

The CardiAct CC-90 Alarmed Outdoor Cabinet provides a compact and secure cabinet for storage of most AEDs. Designed for confined spaces, with its compact size the cabinet also provides clear visibility of the AED with a transparent polycarbonate door that includes AED signage for easy identification in emergencies. Defending against the weather, the CC-90 provides weatherproof protection against dust, rain, cold, rust and humidity with an IP rating of IP54. The CC-90 is suitable for use with the CardiAct Cabinet Stand for versatile mounting options.

The CC-90s Yellow design enhances visibility making it an ideal choice for AEDs located in the public and accessible to anyone. This cabinet ensures that life-saving Automated External Defibrillators (AEDs) are easily accessible, promoting rapid response.

IP Rating: IP54

CARDIACT Alarmed Outdoor AED Cabinet

