



ADDITIONAL INFORMATION

Weight	2 kg
Dimensions	41 × 33 × 19 cm
Carton Qty	1
Pricing Unit	Each
Size	41 x 33 x 19cm
Min Order Qty	1
UNSPSC	42192416
Tax Schedule	Taxable
Item Length	41
Item Width	33
Item Height	19
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The CardiAct CC-90 Alarmed Outdoor Cabinet provides a compact and secure cabinet for storage of most AEDs. Designed for confined spaces, with its compact size the cabinet also provides clear visibility of the AED with a transparent polycarbonate door that includes ANZCOR approved signage for easy identification in emergencies.

Defending against the weather, the CC-90 provides weatherproof protection against dust, rain, cold, rust and humidity with an IP rating of IP54.

The CC-90 is suitable for use with the CardiAct Cabinet Stand for versatile mounting options.

IP Rating: IP54

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



CARDIACT ALARMED AED CABINET 41 X 33 X 19CM

CC-90

CARDIACT Alarmed Outdoor AED Cabinet