



ADDITIONAL INFORMATION

| | |
|-----------------------------|---------------------|
| Weight | 0.015 kg |
| Dimensions | 8.3 × 3.5 × 11.5 cm |
| Carton Qty | 20 |
| Pricing Unit | Each |
| Min Order Qty | 1 |
| ARTG ID | 386216 |
| GMDN Code | 63222 |
| Medical Device Class | Class 1 |
| Tax Schedule | Non-Taxable |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | No |

PRODUCT DESCRIPTION

CPRWrap® Aid is an innovative tool that helps make CPR easier. In an emergency, just remove the product from its quick-open package, lay it on the chest of the sudden cardiac arrest (SCA) victim, and follow the simple instructions. Made with a non-woven fabric that is eco-friendly, this is more than a CPR mask. It provides visual prompts to guide the rescuer through the AHA-recommended steps of CPR. To ensure sanitary protection, the non-woven chest overlay includes a mouth-to-mouth, one-way breathing barrier, and face shield.

- Step-by-step guide for the four AHA-recommended steps
- Visual guides ensure proper hand placement
- One-way valve mouth barrier protects from contaminants
- Non-woven (eco friendly) fabric
- Easy-to-open package for quick access in a medical emergency
- AED compatible

- Bilingual (ENG/SPAN)
- For ages 1 - 8

| Aid Size | Product Dimensions | Pouch Measurements | Pouch Weight |
|----------|--------------------|--------------------|--------------|
| Adult | 16" W x 31.5" H | 3.5" W x 3.5" H | 1.4 oz |
| Child | 9.4" W x 23" H | 3.5" W x 3.5" H | 0.8 oz |
| Infant | 8" W x 14.8" H | 3.5" W x 3.5" H | 0.4 oz |

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.