

CARDIACT WATERPROOF TOUGH AED CASE 33 X 28 X 12CM

CW-10





ADDITIONAL INFORMATION

Weight 1.72 kg

Dimensions $34.4 \times 29.3 \times 13.5 \text{ cm}$

Carton Qty 6

Pricing Unit Unit/1

Size 330mm x 280mm x 120mm

Min Order Qty 1

UNSPSC 42192416
Tax Schedule Taxable

Item Length 33

Item Width 28

Item Height 12

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

This sturdy and rugged water and airtight case is ideal for exceptionally humid or rugged conditions. This customisable case is essential for all marine and outdoor applications.

Rated to IP67 this case can also be equipped with the Desiccant reusable Silica Gel to remove any moisture that may be present in the atmosphere after opening the case.

Fits all HeartSine Defibrillators (PAD-500P, PAD-350P, PAD-360P).

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

AED Case for Heartsine PAD-350P/360P/500P/PDU 399. NOTE: Will NOT FIT LifePak AEDs - please see CW-20.



CARDIACT WATERPROOF TOUGH AED CASE 33 X 28 X 12CM

CW-10