



ADDITIONAL INFORMATION

Weight	0.06 kg
Dimensions	13 × 8 × 1 cm
Carton Qty	1
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The reusable silica gel removes any moisture that might be present in the atmosphere once the defibrillator case is opened.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.