

# **AEROWASH DROP EYEWASH BOTTLE 100ML**

AWDREF100



## ADDITIONAL INFORMATION

Weight	0.16 kg
Dimensions	4 × 4 × 21 cm
Carton Qty	25
Pricing Unit	Unit/1
Size	100ml
Min Order Qty	25
ARTG ID	142720
GMDN Code	37207
Medical Device Class	Class 2a
UNSPSC	46181812
Tax Schedule	Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

#### **PRODUCT DESCRIPTION**

The innovative DROP eyewash station refill bottles all fit into the DROP eyewash station (<u>AWDS01</u>), which is a leader in innovation with a drop down dust cover, integrated mirror, integrated eye cup and a secure safety seal. DROP eyewash refills are a popular solution for convenient eye and wound irrigation.

#### **Strictly Sodium Chloride**

• DROP Eyewash is made up strictly of Sodium Chloride 0.9% sterile solution which is safe and cleansing for eyes and wounds

#### Secure Twist top opening

• DROP Eyewash refills have a secure twist top opening that keeps the saline sterile and secure but makes it easily available when necessary.

#### Instructions

• Instructions for use are clearly indicated on the DROP Eyewash saline bottles



## **AEROWASH DROP EYEWASH BOTTLE 100ML**

AWDREF100

### Also Available in 250ml and 500ml Refills

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Drop Eyewash Refill