

AWDS01



ADDITIONAL INFORMATION

Weight 1.8 kg

Dimensions $33 \times 8 \times 41 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Size 500ml x 2

Min Order Qty 1

ARTG ID 142720 **GMDN Code** 37207

Medical Device Class 2a

UNSPSC 46181810

Tax Schedule Taxable

Is Medical Device Yes
Is Serialized Item No
Is Dangerous Good No
Shelf Life (Months) 36

PRODUCT DESCRIPTION

The innovative DROP eyewash station is a leader in innovation with a drop down dust cover, integrated mirror, integrated eye cup and a secure safety seal.

Easy Refills

DROP eyewash refills are a popular solution for convenient eye and wound irrigation. Available in 500ml, 250ml and 100ml refills.

Strictly Sodium Chloride



AWDS01

	 DROP Eyewash is made up strictly of Sodium Chloride 0.9% sterile solution which is safe and cleansing for eyes and wounds
Wall mour	ntable
	o This station allows an eye wash resource to be mounted so it is easy seen and available for emergency use.
	
Mirror	
	This built in mirror allows a victim to see debris in the eye and self administer first aid
Built in eye bath	
	o DROP Eyewash have a built in eye bath in the top of the bottle for an all in one eye cleansing operation



AWDS01

Secure Twist top opening	
 DROP Eyewash refills have a secure twist top opening that keeps the saline sterile and secure but makes it easily available when necessary. 	
Instructions	
Instructions for use are clearly indicated on the DROP Eyewash saline bottles	
Clear Dust Cover	
o The clear lid on the DROP Eyewash Station keeps it free of dust and allows the contents to be clearly visible.	
This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.	



AWDS01

Eyewash Station