



ADDITIONAL INFORMATION

Weight	1.8 kg
Dimensions	33 × 8 × 41 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	500ml x 2
Min Order Qty	1
ARTG ID	142720
GMDN Code	37207
Medical Device Class	Class 2a
UNSPSC Code	46181810
Tax Schedule	Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36



PRODUCT DESCRIPTION

The AEROWASH Drop Eyewash Station provides an immediately accessible irrigation point for workplaces and public facilities. The wall-mount design presents two ready-to-use 500 mL bottles of sterile 0.9% sodium chloride to assist with clearing and cleaning the ocular area or rinsing minor wounds. A clear, drop-down dust cover protects the contents while keeping visibility high, so responders can locate and access the unit quickly when seconds matter.

Thoughtfully integrated features support calm, accurate first response. The built-in mirror helps users or attendants check the affected eye and align the rinse stream, while the moulded eye cup on each bottle provides a more directed, comfortable flow. A secure safety seal maintains confidence that solution remains sterile until required, and the compact footprint suits first aid rooms, workshops, laboratories, and vehicles.



AEROWASH DROP EYEWASH STATION (2 X 500ML BOTTLES)

AWDS01

Highly visible and easy to install, it is a practical addition to compliance-focused environments that value rapid access to irrigation.

- **Ready when needed:** wall-mount design with clear drop-down dust cover.
- **Adequate irrigation volume:** two 500 mL bottles of sterile 0.9% sodium chloride.
- **Faster assessment:** integrated mirror supports quick alignment of the rinse stream.
- **Directed flow:** built-in eye cup promotes comfortable, focused irrigation.
- **Simple upkeep:** compatible with 100/250/500 mL DROP eyewash refills.

Includes two 500 mL bottles of sterile 0.9% sodium chloride, providing ample irrigation volume for immediate response while keeping replacements simple and predictable for busy workplaces.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Eyewash Station