



ADDITIONAL INFORMATION

Weight	1.8 kg
Dimensions	33 × 8 × 41 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	500ml x 2
Min Order Qty	1
ARTG ID	142720
GMDN Code	37207
Medical Device Class	Class 2a
UNSPSC	46181810
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394004191

PRODUCT DESCRIPTION

The innovative DROP eyewash station is a leader in innovation with a drop down dust cover, integrated mirror, integrated eye cup and a secure safety seal.

Easy Refills

DROP eyewash refills are a popular solution for convenient eye and wound irrigation. Available in [500ml](#), [250ml](#) and [100ml](#) refills.

Strictly Sodium Chloride

- DROP Eyewash is made up strictly of Sodium Chloride 0.9% sterile solution which is safe and cleansing for eyes and wounds

Wall mountable

- This station allows an eye wash resource to be mounted so it is easy seen and available for emergency use.

Mirror

- This built in mirror allows a victim to see debris in the eye and self administer first aid

Built in eye bath

- DROP Eyewash have a built in eye bath in the top of the bottle for an all in one eye cleansing operation

Secure Twist top opening

- DROP Eyewash refills have a secure twist top opening that keeps the saline sterile and secure but makes it easily available when necessary.

Instructions

- Instructions for use are clearly indicated on the DROP Eyewash saline bottles

Clear Dust Cover

- The clear lid on the DROP Eyewash Station keeps it free of dust and allows the contents to be clearly visible.

Eyewash Station



AEROWASH DROP EYEWASH STATION (2 X 500ML BOTTLES)

AWDS01