



ADDITIONAL INFORMATION

Weight	1.8 kg
Dimensions	33 × 8 × 41 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	500ml x 2
Min Order Qty	1
ARTG ID	142720
GMDN Code	37207
Medical Device Class	Class 2a
UNSPSC	46181810
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	36

PRODUCT DESCRIPTION

The innovative DROP eyewash station is a leader in innovation with a drop down dust cover, integrated mirror, integrated eye cup and a secure safety seal.

Easy Refills

DROP eyewash refills are a popular solution for convenient eye and wound irrigation. Available in [500ml](#), [250ml](#) and [100ml](#) refills.

Strictly Sodium Chloride

- DROP Eyewash is made up strictly of Sodium Chloride 0.9% sterile solution which is safe and cleansing for eyes and wounds

Wall mountable

- This station allows an eye wash resource to be mounted so it is easy seen and available for emergency use.

Mirror

- This built in mirror allows a victim to see debris in the eye and self administer first aid

Built in eye bath

- DROP Eyewash have a built in eye bath in the top of the bottle for an all in one eye cleansing operation

Secure Twist top opening

- DROP Eyewash refills have a secure twist top opening that keeps the saline sterile and secure but makes it easily available when necessary.

Instructions

- Instructions for use are clearly indicated on the DROP Eyewash saline bottles

Clear Dust Cover

- The clear lid on the DROP Eyewash Station keeps it free of dust and allows the contents to be clearly visible.

Eyewash Station



AEROWASH DROP EYEWASH STATION (2 X 500ML BOTTLES)

AWDS01