

AERORESCUE ALLOY DUAL-FOLD EMERGENCY POLE STRETCHER

STR-02

ADDITIONAL INFORMATION



Weight 6.35 kg

Dimensions $109 \times 20 \times 10 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Size 2070mm x 540mm x 140mm

Min Order Qty 1

UNSPSC 42192207
Tax Schedule Taxable

Item Length207Item Width54Item Height14Item Measurement UnitcmIs Serialized ItemNo

Is Dangerous Good No

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-02

Weight limit: 160kg

Stretcher Weight: 5.5kg

Popular basic stretcher offering great value for money. Folds to 103cm x 18cm x 10cm in a convenient carry bag

Stretcher material is wipeable (vinyl-like). Carry bag is non-wipeable.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AERORESCUE ALLOY DUAL-FOLD EMERGENCY POLE STRETCHER

STR-02

Alloy Emergency Pole Stretcher