



## ADDITIONAL INFORMATION

<b>Weight</b>	6.75 kg
<b>Dimensions</b>	70 × 32 × 17 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	2140mm x 540mm x 130mm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236801
<b>GMDN Code</b>	13818
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42192207
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	214
<b>Item Width</b>	54
<b>Item Height</b>	13
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

### STR-01

**Weight limit:** 160kg

Folds into a compact carry case measuring 64cm x 27cm x 13cm

Stretcher material is non-wipeable (canvas-like). Case is non-wipeable.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



# AERORESCUE ALLOY QUAD-FOLD EMERGENCY POLE STRETCHER WITH CARRY CASE

STR-01

Alloy Emergency Pole Stretcher