



## ADDITIONAL INFORMATION

<b>Weight</b>	5 kg
<b>Dimensions</b>	213 × 10 × 20 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	2070mm x 540mm x 130mm
<b>Min Order Qty</b>	1
<b>UNSPSC</b>	42192207
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	207
<b>Item Width</b>	54
<b>Item Height</b>	13
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	29341394019519

## PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

### STR-03

**Weight limit:** 160kg

**Stretcher Weight:** 4.5kg

Product Size: 215cm x 55cm x 13cm. Folds down to 182cm x 17cm x 9cm.

Popular single fold stretcher offering great value for money.



# AERORESCUE ALLOY SINGLE-FOLD EMERGENCY POLE STRETCHER (FOLDS WIDTH WAYS)

STR-03

Alloy Emergency Pole Stretcher