



## ADDITIONAL INFORMATION

<b>Weight</b>	8.7 kg
<b>Dimensions</b>	172 × 6 × 47 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	1900mm x 430mm x 60mm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236801
<b>GMDN Code</b>	13818
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42171612
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	190
<b>Item Width</b>	43
<b>Item Height</b>	6
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

### STR-06

**Weight limit:** 160kg

Scoop stretcher designed for gently scooping up a patient while maintaining the found position minimising the risk of further injury. Extendable head rest for 170cm-190cm patient size.

Folds down to 172cm x 43cm x 6cm.

Does NOT include a carry case/bag.

NOTE: The stretcher is supplied in either yellow or silver and the colour provided will vary based on current stock levels.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Alloy Scoop Stretcher

NOTE: The stretcher is supplied in either yellow or silver and the colour provided will vary based on current stock levels.