

ADDITIONAL INFORMATION



Weight	6 kg
Dimensions	99 × 10 × 56 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	1850mm x 515mm x 210mm
Min Order Qty	1
ARTG ID	236801
GMDN Code	13818
Medical Device Class	Class 1
UNSPSC Code	42192207
Tax Schedule	Taxable
Item Length	185
Item Width	51.5
Item Height	21
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-05:

Weight limit: 160kg

Stretcher with wheels and handles for emergency patient transport.

Folds to 93cm x 52cm x 9cm

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AERORESCUE ALLOY FOLDAWAY EMERGENCY STRETCHER WITH WHEELS

STR-05

Alloy Emergency Foldaway Stretcher