

ADDITIONAL INFORMATION

Weight	6 kg
Dimensions	99 × 10 × 56 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	1850mm x 515mm x 210mm
Min Order Qty	1
UNSPSC	42192207
Tax Schedule	Taxable
Item Length	185
Item Width	51.5
Item Height	21
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	29341394019526



PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-05:

Weight limit: 160kg

Stretcher with wheels and handles for emergency patient transport.

Folds to 93cm x 52cm x 9cm

Alloy Emergency Foldaway Stretcher



AERORESCUE ALLOY FOLDAWAY EMERGENCY STRETCHER WITH WHEELS

STR-05