



## ADDITIONAL INFORMATION

<b>Weight</b>	0.06 kg
<b>Dimensions</b>	17.5 × 10.6 × 3.4 cm
<b>Carton Qty</b>	150
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	10cm x 17cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	201258
<b>GMDN Code</b>	201258
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	17
<b>Item Width</b>	10
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

FirstCare Trauma & Haemorrhage Control Bandage consolidates a sterile wound pad, robust elastic wrap, and an integral pressure applicator into a single, ready-to-deploy solution for external bleeding control. Built for pre-hospital responders, workplace first aiders, and community emergencies, it enables fast placement of direct pressure and helps maintain it through movement and handover. The all-in-one format reduces the number of separate pieces to carry and simplifies standard operating procedures when time is critical.

The sterile pad is positioned directly over the site while the wrap provides circumferential tension. A built-in pressure applicator

focuses compression exactly where it is needed to address bleeding, and a secondary dressing component helps keep the pad stable as the patient is moved. The non-adherent contact layer is designed for cleaner release at changeout, supporting low-trauma removal and efficient progression to definitive care.

To streamline application in confined spaces or stressful conditions, the system includes a closure bar that fastens without clips, pins, tape, or knots. The design also supports one-handed self-application, aiding effective use when assistance is limited. This practical construction supports consistent technique across teams, reduces training complexity, and helps responders establish targeted compression quickly with minimal equipment.

- **Direct pressure made simple:** integral applicator concentrates compression at the wound.
- **Sterile pad with secondary dressing:** helps maintain placement and pressure in transit.
- **Fast, tool-free closure:** built-in bar avoids clips, pins, tape, and knots.
- **Low-trauma removal:** non-adherent contact layer designed to lift cleanly.
- **Supports one-handed use:** aids effective application when help is limited.

The **Size: 10 x 17 cm** option provides a focused sterile pad area with ample wrap to establish targeted compression, balancing coverage and manoeuvrability for rapid bleeding control.

Finished in **Military Green** to suit low-profile field kits, this configuration pairs the 10 x 17 cm pad with the integral pressure applicator for focused compression while maintaining a subdued, professional appearance.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Military Trauma & Haemorrhage Control Wound Dressing (green)