

FIRSTCARE CIVILIAN TRAUMA & HEMORRHAGE CONTROL BANDAGE 10 X 17CM (WHITE)

FCP-03



ADDITIONAL INFORMATION

Weight 0.06 kg

Dimensions $15 \times 3 \times 6.5 \text{ cm}$

Carton Qty 150

Pricing Unit Unit/1

Size 10cm x 17cm

Min Order Qty 1

ARTG ID 201258
GMDN Code 201258
Medical Device Class Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Item Length 17

Item Width 10

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

FirstCare Trauma & Hemorrhage Control Bandage streamlines emergency bleeding control by combining a sterile wound pad, a robust wrap, and an integral pressure applicator into one compact solution. Designed for pre-hospital responders and workplace first aid, it helps apply direct pressure rapidly and keep it in place during patient movement and handover.

The sterile pad is positioned directly over the wound, while the wrap secures tension across the site. A built-in pressure applicator concentrates force where it is needed to address bleeding. Once positioned, a secondary dressing component maintains pad placement and pressure, supporting efficient care until advanced treatment is available.

To simplify application when time is critical, the bandage features a closure bar that fastens without clips, pins, tape, or knots. The non-adherent pad is designed to release cleanly at changeout, and one-handed self-application is supported for scenarios where assistance may be limited. This practical, all-in-one configuration reduces the number of separate items to manage and train on, helping standardise response across teams.



FIRSTCARE CIVILIAN TRAUMA & HEMORRHAGE CONTROL BANDAGE 10 X 17CM (WHITE)

FCP-03

- **Direct pressure made simple:** integral applicator focuses compression at the wound.
- Sterile pad with secondary dressing: helps secure pressure and placement in transit.
- Fast, tool-free closure: built-in bar eliminates pins, clips, tape, and knots.
- **Low-trauma removal:** non-adherent pad designed to lift cleanly.
- **Supports one-handed use:** aids effective application when help is limited.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Civilian Trauma & Haemorrhage Control Wound Dressing (white)