

**ADDITIONAL INFORMATION**

<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	257365
<b>GMDN Code</b>	42424
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42271913
<b>Tax Schedule</b>	Non-Taxable
<b>Location</b>	AU-Cameron Dr
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



**PRODUCT DESCRIPTION**

The Guedel Airway set of 6 is an essential tool for maintaining an open airway in unconscious patients. The Guedel Airway prevents occlusion due to tongue relaxation and feature a central lumen that accommodates various suction catheters. The airways are color-coded for easy size identification and feature a thickened area that acts as a bite block, protecting the airway from biting. The set includes airways in sizes ranging from neonatal to large adult, making it a versatile and comprehensive solution for a variety of patients.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Guedel Airway Set of 6 (#0 - #5)