

ADDITIONAL INFORMATION

Pricing Unit	Each
Min Order Qty	1
ARTG ID	257365
GMDN Code	42424
Medical Device Class	Class 1s
UNSPSC	42271913
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The Guedel Airway set of 6 is an essential tool for maintaining an open airway in unconscious patients. The Guedel Airway prevents occlusion due to tongue relaxation and feature a central lumen that accommodates various suction catheters. The airways are color-coded for easy size identification and feature a thickened area that acts as a bite block, protecting the airway from biting. The set includes airways in sizes ranging from neonatal to large adult, making it a versatile and comprehensive solution for a variety of patients.

Guedel Airway Set of 6 (#0 - #5)