

GUEDEL AIRWAY SIZE 0 - SIZE 5 PACK/6

GAB-SET



ADDITIONAL INFORMATION

Pricing Unit Each
Min Order Qty 1

ARTG ID 257365

GMDN Code 42424

42424

Medical Device Class

UNSPSC 42271913

Tax Schedule Non-Taxable

Location AU-Cameron Dr

Class 1s

Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

The Guedel Airway set of 6 is an essential tool for maintaining an open airway in unconscious patients. The Guedel Airway prevents occlusion due to tongue relaxation and feature a central lumen that accommodates various suction catheters. The airways are color-coded for easy size identification and feature a thickened area that acts as a bite block, protecting the airway from biting. The set includes airways in sizes ranging from neonatal to large adult, making it a versatile and comprehensive solution for a variety of patients.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Guedel Airway Set of 6 (#0 - #5)