



ADDITIONAL INFORMATION

Weight	0.325 kg
Dimensions	19.5 × 13 × 3.5 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	Paediatric
Min Order Qty	1
ARTG ID	312730
GMDN Code	47912
Medical Device Class	Class 2b
UNSPSC Code	42172105
Tax Schedule	Non-Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	Yes
Shelf Life (Months)	48
Warranty Period (Months)	12

PRODUCT DESCRIPTION

The HeartSine Pad-Paks are designed to provide a convenient and reliable solution for users of the HeartSine defibrillator series. These innovative cartridges contain both electrodes and batteries in one cassette, eliminating the need to worry about expiring components. This concept is not only more efficient but also more cost-effective throughout the device's lifecycle. The Pad-Paks are suitable for use with the HeartSine defibrillator series, including the PAD 350P, PAD 360P, and PAD 500.

The cartridges have a shelf life of four years from the date of manufacture and provide 60 shocks or continuous monitoring over a period of six hours. Installation is straightforward, requiring just one click to install the cartridge. This simplicity and ease of use make the HeartSine Pad-Paks an ideal choice for those who require a reliable and efficient defibrillator solution.

The HeartSine Pad-Paks are a single use device, containing a battery power source and defibrillation electrodes that can only be used with the samaritan PAD devices.

The Paediatric size is designed for use on children under 25kg.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



HEARTSINE PINK PAD-PAK PADS & BATTERY PACK - PAEDIATRIC (DG)

PAD-PAK-04

directions for use.

PAD-Pak (Pink)