



ADDITIONAL INFORMATION

Weight	0.22 kg
Dimensions	29.5 × 9.9 × 4.8 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	Set of 25
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
Tax Schedule	Taxable
Item Length	4.8
Item Width	9.9
Item Height	29.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The easy-to-use samaritan PAD Trainer with remote control guides users through simulated analysis, energy delivery and prompted Cardio Pulmonary Resuscitation (CPR) intervals with the look and feel of a live samaritan Public Access Defibrillator (PAD) without the actual charge and discharge of an electrical shock. Available in 500P, 350P and 360P configurations, the samaritan training device is designed to help and support trainers fulfil their training obligations.

HeartSine encourage everyone to take part in a training course to help combat the loss of life through SCA. The Universal Trainer is one device designed to help / support the trainers fulfil their training obligations.

[AED Website](#)



HEARTSINE TRAINER DEFIBRILLATOR PADS PACK/25

TRN-ACC-03