



ADDITIONAL INFORMATION

Weight	0.1 kg
Dimensions	15 × 11.5 × 3 cm
Pricing Unit	Unit/1
Size	One Size
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	JC10105060167124495

PRODUCT DESCRIPTION

The easy-to-use samaritan PAD Trainer with remote control guides users through simulated analysis, energy delivery and prompted Cardio Pulmonary Resuscitation (CPR) intervals with the look and feel of a live samaritan Public Access Defibrillator (PAD) without the actual charge and discharge of an electrical shock. Available in 500P, 350P and 360P configurations, the samaritan training device is designed to help and support trainers fulfil their training obligations.

HeartSine encourage everyone to take part in a training course to help combat the loss of life through SCA. The Universal Trainer is one device designed to help / support the trainers fulfil their training obligations.

[AED Website](#)

Trainer Pad-Pak Samaritan wind-up