



ADDITIONAL INFORMATION

Weight	1.5 kg
Dimensions	29 × 23.5 × 14.8 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	One Size
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The easy-to-use samaritan PAD Trainer with remote control guides users through simulated analysis, energy delivery and prompted Cardio Pulmonary Resuscitation (CPR) intervals with the look and feel of a live samaritan Public Access Defibrillator (PAD) without the actual charge and discharge of an electrical shock. Available in 500P, 350P and 360P configurations, the samaritan training device is designed to help and support trainers fulfil their training obligations.

HeartSine encourage everyone to take part in a training course to help combat the loss of life through SCA. The Universal Trainer is one device designed to help / support the trainers fulfil their training obligations.

[AED Website](#)

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



HEARTSINE SAMARITAN 500P TRAINER DEFIBRILLATOR

TRN-500P

Defibrillator Trainer Samaritan PAD 500P