



ADDITIONAL INFORMATION

Weight	0.36 kg
Dimensions	21.7 × 4 × 13.9 cm
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Non-Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

- Child / Infant pads for SCA victims under 8 y.o., or weighing less than 25kg.
- Packaged in a cartridge that can clips into the AED
- Adhesive pads
- Clearly printed diagrams to show rescuers the correct pad placement, as well as sensors to indicate when the pads have been place correctly, improving the effectiveness of the defibrillator's shock.
- 2 year shelf life
- Compatible with M5066A

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



HEARTSTART HS1 CHILD SMART PADS PAIR/2

HS1-PAD-PED