



ADDITIONAL INFORMATION

Weight	0.28125 kg
Dimensions	4.2 × 10.5 × 11 cm
Carton Qty	32
Pricing Unit	Box/1
Size	23.5cm x 12cm
Min Order Qty	32
ARTG ID	159730
GMDN Code	34964
Medical Device Class	Class 2a
UNSPSC Code	42142108
Tax Schedule	Taxable
Item Length	23.5
Item Width	12
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AeroCool Instant Ice Pack is a portable, single-use cold therapy solution designed to deliver rapid cooling wherever it is needed. The compact pouch is ready to activate on the spot, making it a practical choice for sports sidelines, workplaces, schools, clinics and home first aid kits. Once activated, the pack provides a soothing cold sensation that can be applied comfortably to the body surface.

The formulation uses urea and water to generate instant cooling without the need for pre-chilling or refrigeration. This makes AeroCool Instant Ice Pack particularly useful in environments where access to freezers is limited or when unexpected injuries occur

away from fixed facilities. The flexible pack conforms to different body areas, supporting contact around joints and contoured surfaces.

AeroCool Instant Ice Pack is suitable for a wide range of everyday bumps and strains. It is used to help reduce localised discomfort and swelling from minor sprains, muscular strain, bruising, knocks and bumps. It can also be used as part of general first aid management where temporary cooling is recommended, such as helping to lower temperature from fevers, easing toothache discomfort or soothing itchiness from insect bites.

The lightweight construction and streamlined pouch profile allow the pack to fit easily into first aid kits, sports bags, glove boxes and desk drawers. With simple activation and no requirement for power or pre-cooling, AeroCool Instant Ice Pack supports quick response in busy or remote settings and is a practical addition to hot and cold treatment ranges.

- **Instant cooling:** chemical reaction delivers cold therapy without refrigeration.
- **Portable format:** compact pouch is easy to store in kits, bags and vehicles.
- **Versatile applications:** suitable for minor sprains, strains, bruises and everyday knocks.
- **Comfortable contact:** flexible pack conforms to different body areas.
- **Ready when needed:** ideal for environments where traditional ice is not available.

Available in both 80g and 240g sizes. This larger 240g size is great for bruises and larger knocks and bumps.

This pack is supplied as a single ice pack in a box, offering a convenient individual unit that fits easily into first aid kits, sports bags or workplace cabinets for quick access when cooling is needed.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Instant Ice Pack 240g