

**ADDITIONAL INFORMATION**

Weight	7 kg
Carton Qty	1
Pricing Unit	Unit/1
Size	110 x 50 x 90.5
Min Order Qty	1
UNSPSC Code	24101643
Tax Schedule	Taxable
Item Length	110
Item Width	50
Item Height	90.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION**Medical/Instrument Trolleys**

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- 2 shelves
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking



LARGE STAINLESS STEEL TROLLEY 110 X 50 X 90.5CM

ACR10

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Large Trolley