



## ADDITIONAL INFORMATION

<b>Weight</b>	2 kg
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Non-Taxable
<b>Location</b>	AU-Cameron Dr
<b>Is Serialized Item</b>	Yes
<b>Is Dangerous Good</b>	Yes

## PRODUCT DESCRIPTION

- Rigorously drop-tested device
- The LIFEPAK is programmable, with adjustable settings to match your team's CPR and resuscitation protocols
- The AED is easy-to-use, with loud voice prompts and on-screen graphics providing guidance on applying electrodes and initiating a shock - Vehicle friendly designed to ride along without damage from continuous vibrations and other movements
- Flexible power options between a rechargeable Li-ion or non-rechargeable Lithium-Manganese Dioxide battery
- Escalates defibrillation power to an industry-leading 360J
- CPRMAX™ Technology minimized CPR interruptions by allowing compressions to continue during AED charging
- CPR Countdown Timer provides direction for the length of hands-on time for each CPR period based on system protocol
- 3-lead ECG function is available when need to provide accurate heart rhythm monitoring
- Shock counter digitally records and displays delivered shocks for added insight
- Large, easy-to-read LCD screen provides more information at a glance
- Electrodes a fully compatible with all other LIFEPAK defibrillators and monitors

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.