



ADDITIONAL INFORMATION

Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The Trainer 1000 provides realistic training for LIFEPAK 1000 users. It utilizes the same screen messages, audible tones and voice prompts as those found in the LIFEPAK 1000 Defibrillator to guide users-in-training through simulated analysis, energy delivery and prompted CPR intervals.

- Mimics the LIFEPAK 1000 defibrillator in appearance and operation
- Simulates AED and Manual Mode
- Includes simulated cprMAX™ technology
- Wireless remote control enables flexible training scenarios

A versatile trainer, the Trainer 1000 simulates both AED and Manual Mode. All functions are programmable to simulate the full range of possible set-up options found in the LIFEPAK 1000 Defibrillator. The Trainer 1000 allows you to run preconfigured shock scenarios or create a dynamic custom training event of your own using the included wireless remote control. With a rechargeable Lithium-ion battery, the Trainer 1000 operates up to 4 hours. Realistic, versatile and easy to use, the Trainer 1000 is ideal for training large groups and as an alternative to using a live defibrillator for training. Mimics the LIFEPAK 1000 defibrillator in appearance and operation Simulates AED and Manual Mode Includes simulated cprMAX™ technology Wireless remote control enables flexible training scenarios.