



ADDITIONAL INFORMATION

| | |
|-----------------------------|--------------------|
| Weight | 0.32 kg |
| Dimensions | 15.5 × 5 × 24.4 cm |
| Carton Qty | 1 |
| Pricing Unit | Each |
| Min Order Qty | 1 |
| ARTG ID | 120810 |
| GMDN Code | 45806 |
| Medical Device Class | Class 1 |
| Tax Schedule | Non-Taxable |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | Yes |
| Shelf Life (Months) | 24 |

PRODUCT DESCRIPTION

LIFEPAK CR Plus Charge Pak contains the Pads and Battery Pak for the LifePak CR-Plus AED and the LifePak Express AED.

Includes 1 x Pair of Adult Pads & 1 x Battery Pak. The Charge Pak has a two-year expiry.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.