



## ADDITIONAL INFORMATION

<b>Weight</b>	0.22 kg
<b>Dimensions</b>	16.5 × 3.9 × 26.7 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	232970
<b>GMDN Code</b>	41587
<b>Medical Device Class</b>	Class 2b
<b>Tax Schedule</b>	Non-Taxable
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	216

## PRODUCT DESCRIPTION

The LIFEPAK CR Plus Paediatric Pads are a must-have for any first responder or medical professional who may need to treat infants or children in emergency situations. These pads are designed specifically for the Lifepak CR Plus and Lifepak 1000 AED monitors and are a replacement for the standard adult pads.

The paediatric defibrillation pads are designed to provide reduced energy defibrillation, making them safe and effective for use on infants and children under the age of 8 or weighing less than 25kg. This ensures that the treatment is tailored to the unique needs of paediatric patients, minimizing the risk of injury or harm.

Using the paediatric pads on adult patient, will not have any negative effects on the adult but it will not be efficient for the treatment of adult.

These replacement pads are easy to install, and their compact design makes them easy to transport and store. Whether you're a first responder, a medical professional, or a parent, the LIFEPAK CR Plus Paediatric Pads are an essential piece of equipment for ensuring the safety and well-being of children in emergency situations.

The pads have an 18-month life.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

