

# LIFEPAK CR2 ESSENTIAL SEMI-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E



### **ADDITIONAL INFORMATION**

Weight 4.02 kg

**Dimensions**  $31.5 \times 25.5 \times 19 \text{ cm}$ 

Carton Qty 1

Pricing Unit Unit/1

**Size** One Size

Min Order Qty 1

ARTG ID 282772

Medical Device Class 2b

UNSPSC 42172101

Tax Schedule Non-Taxable

Is Serialized Item Yes
Is Dangerous Good Yes

#### **PRODUCT DESCRIPTION**

Sudden cardiac arrest (SCA) can happen to anyone, anywhere. Immediate treatment is vital. A victim 's chance of survival dramatically decreases for every minute without treatment. That's why public access defibrillators are so important: they put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best.

## Layered design

Layered design with easy to follow bold graphics. Both trained and untrained AED users clearly know how to begin.

### **QUIK-STEP**<sup>™</sup> electrodes



Semi-Automatic Defibrillator AED

# LIFEPAK CR2 ESSENTIAL SEMI-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E

Peel directly off the base for faster placement.
Metronome and CPR coaching
Sets an effective pace and audibly guides users.
ClearVoice <sup>™</sup> technology
Detects background noise and adjusts volume accordingly for clear instruction.
Fully automatic
Available in fully or semi-automatic models.
Highest available escalating energy
Up to 360J for more effective shocks as needed.
Bilingual
Toggle between two pre-set languages when using the device.
Child mode
Child mode delivers lower energy levels appropriate for young children without having to change electrodes.
LIFEPAK TOUGH™
IP55 rating for challenging environments.
Battery Expiry (CR-2-BAT): 4-year from install date, with 1-year installation window from date of manufacture.
<b>Electrode Expiry:</b> 4-years from manufacutring date. This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



# LIFEPAK CR2 ESSENTIAL SEMI-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E