



ADDITIONAL INFORMATION

Weight	4.02 kg
Dimensions	31.5 × 25.5 × 19 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	One Size
Min Order Qty	1
ARTG ID	282772
Medical Device Class	Class 2b
UNSPSC Code	42172101
Tax Schedule	Non-Taxable
Location	AU-Cameron Dr
Is Medical Device	Yes
Is Serialized Item	Yes
Is Dangerous Good	Yes

PRODUCT DESCRIPTION

This semi-automatic configuration places control of shock delivery in the responder's hands, combining clear voice prompts and a flashing shock button so users can visually confirm that everyone is clear before activating therapy.

Supplied as a single defibrillator unit, this LIFEPAK CR2 Essential Semi-Automatic model is ideal where one dependable device is needed for a workplace, community venue, or household, combining advanced analysis technology and CPR coaching in a compact, easy-to-recognise design.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Semi-Automatic Defibrillator AED



LIFEPAK CR2 ESSENTIAL SEMI-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E