

LIFEPAK CR2 ESSENTIAL FULLY-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E-A



ADDITIONAL INFORMATION

Weight 4.06 kg

Dimensions $31.5 \times 25.5 \times 19 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Size One Size

Min Order Qty 1

ARTG ID 282772 **GMDN Code** 48047

Medical Device Class 2b

UNSPSC 42172101

Tax Schedule Non-Taxable

Is Serialized Item Yes
Is Dangerous Good Yes

PRODUCT DESCRIPTION

Sudden cardiac arrest (SCA) can happen to anyone, anywhere. Immediate treatment is vital. A victim 's chance of survival dramatically decreases for every minute without treatment. That's why public access defibrillators are so important: they put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best.

Layered design

Layered design with easy to follow bold graphics. Both trained and untrained AED users clearly know how to begin.

QUIK-STEP[™] electrodes

Peel directly off the base for faster placement.

Metronome and CPR coaching

Sets an effective pace and audibly guides users.

ClearVoice[™] technology



LIFEPAK CR2 ESSENTIAL FULLY-AUTOMATIC DEFIBRILLATOR (DG)

CR-2E-A

Detects background noise and adjusts volume accordingly for clear instruction.

Fully automatic

Available in fully or semi-automatic models.

Highest available escalating energy

Up to 360J for more effective shocks as needed.

Bilingual

Toggle between two pre-set languages when using the device.

Child mode

Child mode delivers lower energy levels appropriate for young children without having to change electrodes.

LIFEPAK TOUGH™

IP55 rating for challenging environments.

Battery Expiry (CR-2-BAT): 4-year from install date, with 1-year installation window from date of manufacture.

Electrode Expiry: 4-years from manufacutring date.

Fully-Automatic Defibrillator AED