



## ADDITIONAL INFORMATION

<b>Weight</b>	4.1 kg
<b>Dimensions</b>	31.5 × 25.5 × 19 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	One Size
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	282772
<b>GMDN Code</b>	48047
<b>Medical Device Class</b>	Class 2b
<b>UNSPSC Code</b>	42172101
<b>Tax Schedule</b>	Non-Taxable
<b>Location</b>	AU-Cameron Dr
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	Yes
<b>Is Dangerous Good</b>	Yes

## PRODUCT DESCRIPTION

The LifePak CR2 is designed for simplicity and speed. It's easy to use, even for minimally trained users, thanks to its two-step operation and no-button functionality. Just open the lid, and the device will automatically start, analysing the patient's heart rhythm and delivering a shock if needed—no need to press any buttons.

### Key Features

- **Layered design:** Easy-to-follow bold graphics ensure both trained and untrained users clearly understand how to begin.
- **QUIK-STEP™ electrodes:** Peel directly off the base for faster placement.
- **CPR Insight analysis technology:** Analyses for shockable rhythm during chest compressions without needing to pause.
- **Metronome and CPR coaching:** Sets an effective pace and audibly guides users.
- **ClearVoice™ technology:** Detects background noise and adjusts volume accordingly for clear instruction.
- **Fully automatic or semi-automatic:** Available in fully or semi-automatic models.
- **Highest available escalating energy:** Delivers shocks of up to 360J as needed.
- **Child mode:** Delivers lower energy levels appropriate for young children without the need to change electrodes.
- **LIFEPAK TOUGH™:** IP55 rating makes it suitable for challenging environments.
- **8-year warranty:** Backed by an 8-year warranty for peace of mind.



# LIFEPAK CR2 CPRINSIGHT FULLY-AUTOMATIC DEFIBRILLATOR WITH WI-FI

CR-2-A

## Connectivity

- **Wireless connectivity:** Enables real-time device status monitoring and remote management through the LIFELINKcentral™ AED Program Manager.
- **Automatic notifications:** Immediate alerts when AED is used or requires servicing.
- **Accessory tracking:** Early notifications on expiration of pads and battery to ensure on-time replacement
- **Remote software updates:** Ensures the device maintains optimal functionality through automatic updates.

## Maintenance and Accessories

- **Battery life:** 4-year standby battery life.
  - Battery Expiry (CR-2-BAT): 4 years from install date, with a 1-year installation window from manufacture date.
- **Electrode expiry:** 4 years from manufacturing date.
- **Hard-shell carry case:** Protects the device and enhances portability

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Fully-Automatic Defibrillator AED WiFi