

## LIFEPAK CR2 SEMI-AUTOMATIC DEFIBRILLATOR WITH WI-FI (DG)

CR-2



#### **ADDITIONAL INFORMATION**

Weight 4.04 kg

**Dimensions**  $31.5 \times 25.5 \times 19 \text{ cm}$ 

Carton Qty 1

Pricing Unit Unit/1

**Size** One Size

Min Order Qty 1

ARTG ID 282772

Medical Device Class Class 2b

UNSPSC 42172101

**Tax Schedule** Non-Taxable

Is Serialized Item Yes
Is Dangerous Good Yes

#### PRODUCT DESCRIPTION

The LifePak CR2 is designed for simplicity and speed. It's easy to use, even for minimally trained users, thanks to its two-step operation and no-button functionality. Just open the lid, and the device will automatically start, analysing the patient's heart rhythm and delivering a shock if needed—no need to press any buttons.

#### **Key Features**

- Layered design: Easy-to-follow bold graphics ensure both trained and untrained users clearly understand how to begin.
- **QUIK-STEP™** electrodes: Peel directly off the base for faster placement.
- CPR Insight analysis technology: Analyses for shockable rhythm during chest compressions without needing to pause.
- Metronome and CPR coaching: Sets an effective pace and audibly guides users.
- ClearVoice™ technology: Detects background noise and adjusts volume accordingly for clear instruction.
- Fully automatic or semi-automatic: Available in fully or semi-automatic models.
- **Highest available escalating energy:** Delivers shocks of up to 360J as needed.
- Child mode: Delivers lower energy levels appropriate for young children without the need to change electrodes.
- **LIFEPAK TOUGH™:** IP55 rating makes it suitable for challenging environments.
- **8-year warranty:** Backed by an 8-year warranty for peace of mind.

#### Connectivity

• Wireless connectivity: Enables real-time device status monitoring and remote management through the LIFELINKcentral™



# LIFEPAK CR2 SEMI-AUTOMATIC DEFIBRILLATOR WITH WI-FI (DG)

CR-2

### AED Program Manager.

- Automatic notifications: Immediate alerts when AED is used or requires servicing.
- Accessory tracking: Early notifications on expiration of pads and battery to ensure on-time replacement
- Remote software updates: Ensures the device maintains optimal functionality through automatic updates.

#### Maintenance and Accessories

- Battery life: 4-year standby battery life.
  - Battery Expiry (CR-2-BAT): 4 years from install date, with a 1-year installation window from manufacture date.
- **Electrode expiry:** 4 years from manufacturing date.
- Hard-shell carry case: Protects the device and enhances portability

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Semi-Automatic Defibrillator AED WiFi