



ADDITIONAL INFORMATION

Weight	4.04 kg
Dimensions	31.5 × 25.5 × 19 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	One Size
Min Order Qty	1
ARTG ID	282772
Medical Device Class	Class 2b
UNSPSC Code	42172101
Tax Schedule	Non-Taxable
Location	AU-Cameron Dr
Is Medical Device	Yes
Is Serialized Item	Yes
Is Dangerous Good	Yes

PRODUCT DESCRIPTION

The LifePak CR2 is designed for simplicity and speed. It's easy to use, even for minimally trained users, thanks to its two-step operation and no-button functionality. Just open the lid, and the device will automatically start, analysing the patient's heart rhythm and delivering a shock if needed—no need to press any buttons.

Key Features

- **Layered design:** Easy-to-follow bold graphics ensure both trained and untrained users clearly understand how to begin.
- **QUIK-STEP™ electrodes:** Peel directly off the base for faster placement.
- **CPR Insight analysis technology:** Analyses for shockable rhythm during chest compressions without needing to pause.
- **Metronome and CPR coaching:** Sets an effective pace and audibly guides users.
- **ClearVoice™ technology:** Detects background noise and adjusts volume accordingly for clear instruction.
- **Fully automatic or semi-automatic:** Available in fully or semi-automatic models.
- **Highest available escalating energy:** Delivers shocks of up to 360J as needed.
- **Child mode:** Delivers lower energy levels appropriate for young children without the need to change electrodes.
- **LIFEPAK TOUGH™:** IP55 rating makes it suitable for challenging environments.
- **8-year warranty:** Backed by an 8-year warranty for peace of mind.

Connectivity

- **Wireless connectivity:** Enables real-time device status monitoring and remote management through the LIFELINKcentral™ AED Program Manager.
- **Automatic notifications:** Immediate alerts when AED is used or requires servicing.
- **Accessory tracking:** Early notifications on expiration of pads and battery to ensure on-time replacement
- **Remote software updates:** Ensures the device maintains optimal functionality through automatic updates.

Maintenance and Accessories

- **Battery life:** 4-year standby battery life.
 - Battery Expiry (CR-2-BAT): 4 years from install date, with a 1-year installation window from manufacture date.
- **Electrode expiry:** 4 years from manufacturing date.
- **Hard-shell carry case:** Protects the device and enhances portability

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Semi-Automatic Defibrillator AED WiFi