

# LIFEPAK QUIK-COMBO REDI-PAK ELECTRODES - ADULT (FOR LP1000, LP12, LP15 & LP20E)

## 11996-000017

LP-PAD



### ADDITIONAL INFORMATION

<b>Weight</b>	0.12 kg
<b>Dimensions</b>	25.5 × 26 × 0.5 cm
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	270015
<b>GMDN Code</b>	35035
<b>Medical Device Class</b>	Class 1
<b>Tax Schedule</b>	Non-Taxable
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	216

### PRODUCT DESCRIPTION

In a cardiac emergency, every second counts, which is why these adult electrode pads have been designed to be preconnected prior to use, which eliminates the necessity to change electrode pads during a cardiac emergency, saving time.

With Edge System(TM) technology, even current distribution is achieved, which minimizes the risk of burns to the skin. The single set of electrodes is able to monitor the ECG, provide noninvasive cardioversion, and perform noninvasive pacing and defibrillation.

These disposable QUIK-COMBO electrodes are designed to be used with QUIK-COMBO therapy cables designed specifically for your device, and are intended to be used with a single patient. The lead wire is 2 feet (0.6 m) long. Due to the latex-free electrodes, there is n

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.