



ADDITIONAL INFORMATION

Weight	2.14 kg
Dimensions	54 × 26 × 10 cm
Carton Qty	2
Pricing Unit	Unit/1
Size	120mm diameter Extension length: 1150mm
Min Order Qty	1
UNSPSC Code	42182601
Tax Schedule	Taxable
Item Length	115
Item Width	12
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

Magnifying Lamp

A magnifying lamp is integral to accurately diagnose and treat small injuries and provide sufficient light for casualty treatment. The AeroSupplies LED Magnifying Lamp comes with a table clamp and has a 1200mm reach/extension. It has lens 6 diopter magnification (2.5x) and is has 10 watt LED lights (6000K).

You can purchase a mobile base ([ACR05-MB](#)) or wall mount ([ACR05-WM](#)) separately that are both compatible with this magnifying lamp for different applications.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



LED MAGNIFYING LAMP WITH TABLE CLAMP (12CM DIAMETER, 115CM EXTENSION)

ACR05-ML

Magnifying Lamp with Table Clamp