

## **MANUAL SUCTION PUMP**

AMSP01



## **ADDITIONAL INFORMATION**

Weight 0.38 kg

**Dimensions**  $20.5 \times 17.5 \times 7.8 \text{ cm}$ 

**Pricing Unit** Each

Min Order Qty 1

**UNSPSC** 42144402

Tax Schedule Non-Taxable

Is Medical DeviceNoIs Serialized ItemNoIs Dangerous GoodNo

## PRODUCT DESCRIPTION

Light and portable suction unit. Can be used with one-handed operation which leaves the other hand free to perform important tasks.

Features a knob to adjust the suction pressure (50% can produce 300mmHG vacuum, while 100% can reach 600mmHG).

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Manual Suction Pump