



## ADDITIONAL INFORMATION

<b>Weight</b>	0.38 kg
<b>Dimensions</b>	20.5 × 17.5 × 7.8 cm
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>UNSPSC</b>	42144402
<b>Tax Schedule</b>	Non-Taxable
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

Light and portable suction unit. Can be used with one-handed operation which leaves the other hand free to perform important tasks. Features a knob to adjust the suction pressure (50% can produce 300mmHG vacuum, while 100% can reach 600mmHG).

Manual Suction Pump