



ADDITIONAL INFORMATION

Weight	0.38 kg
Dimensions	20.5 × 17.5 × 7.8 cm
Pricing Unit	Each
Min Order Qty	1
UNSPSC Code	42144402
Tax Schedule	Non-Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Light and portable suction unit. Can be used with one-handed operation which leaves the other hand free to perform important tasks. Features a knob to adjust the suction pressure (50% can produce 300mmHG vacuum, while 100% can reach 600mmHG).

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Manual Suction Pump