



ADDITIONAL INFORMATION

Weight	16 kg
Carton Qty	1
Pricing Unit	Unit/1
Size	60 x 50 x 97
Min Order Qty	1
UNSPSC	24101643
Tax Schedule	Taxable
Item Length	60
Item Width	50
Item Height	97
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- 1 Drawer: 125mm deep
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking
- Full extension telescopic ball bearing drawer runners

Medium Trolley with Drawer