

MEDIUM STAINLESS STEEL TROLLEY WITH DRAWER 60 X 50 X 97CM

ACR13



ADDITIONAL INFORMATION

Weight 16 kg
Carton Qty 1

Pricing Unit Unit/1

Size 60 x 50 x 97

Min Order Qty 1

UNSPSC 24101643

Tax Schedule Taxable

Item Length 60

Item Width 50

Item Height 97

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- 1 Drawer: 125mm deep
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking
- Full extension telescopic ball bearing drawer runners



MEDIUM STAINLESS STEEL TROLLEY WITH DRAWER 60 X 50 X 97CM

ACR13

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Medium Trolley with Drawer