



### ADDITIONAL INFORMATION

<b>Weight</b>	7 kg
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Each
<b>Size</b>	80 x 50 x 97
<b>Min Order Qty</b>	1
<b>UNSPSC</b>	24101643
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	80
<b>Item Width</b>	50
<b>Item Height</b>	97
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

### PRODUCT DESCRIPTION

#### Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

#### Product Features:

- Guard rails to 3 sides of 2 shelves
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking



# MEDIUM STAINLESS STEEL TROLLEY WITH RAILS

## 80 X 50 X 97CM

ACR09

Medium Trolley with Rails