

**ADDITIONAL INFORMATION**

Weight	7 kg
Carton Qty	1
Pricing Unit	Each
Size	80 x 50 x 97
Min Order Qty	1
UNSPSC Code	24101643
Tax Schedule	Taxable
Item Length	80
Item Width	50
Item Height	97
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION**Medical/Instrument Trolleys**

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- Guard rails to 3 sides of 2 shelves
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking



MEDIUM STAINLESS STEEL TROLLEY WITH RAILS

80 X 50 X 97CM

ACR09

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Medium Trolley with Rails