

**ADDITIONAL INFORMATION**

| | |
|------------------------------|----------------|
| Weight | 7 kg |
| Carton Qty | 1 |
| Pricing Unit | Each |
| Size | 80 x 50 x 90.5 |
| Min Order Qty | 1 |
| UNSPSC | 24101643 |
| Tax Schedule | Taxable |
| Item Length | 80 |
| Item Width | 50 |
| Item Height | 90.5 |
| Item Measurement Unit | cm |
| Is Serialized Item | No |
| Is Dangerous Good | No |

PRODUCT DESCRIPTION**Medical/Instrument Trolleys**

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- High quality 304 grade stainless steel
- Safety fold to all edges



MEDIUM STAINLESS STEEL TROLLEY 80 X 50 X 90.5CM

ACR08

Medium Trolley