

## 3 PANEL MOBILE PRIVACY SCREEN 219 X 180 X 15CM

ACR04



## **ADDITIONAL INFORMATION**

Weight 7.5 kg

**Dimensions**  $183 \times 70 \times 10 \text{ cm}$ 

Carton Qty 1

Pricing Unit Unit/1

**Size** 219 x 180 x 15 cm

Min Order Qty 1

UNSPSC 42191607
Tax Schedule Taxable

**Item Length** 219

**Item Width** 180

**Item Height** 15

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

## **PRODUCT DESCRIPTION**

A mobile privacy screen is an essential piece of equipment to ensure treatment of casualties is discreet and professional. Supplied flat packed with minimal assembly required

- Main frame made of steel tube, Pre-treated and epoxy
- Three panels, easy for folding and using
- Curtain cloth (fabric-cotton) can be washed separately
- Six 1.5-inch castors, two with brakes
- Floor clearance is 135mm

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Mobile Privacy Screen



## 3 PANEL MOBILE PRIVACY SCREEN 219 X 180 X 15CM

ACR04